

## 16S Library PCR Protocol with PNA Blockers

*This protocol is a modified version of the Earth Microbiome Project's PCR protocol and the JGI's PNA PCR Protocol.*

**Notes:** There are 25 each of forward and reverse primer barcodes that can be combined to make 625 total unique barcode combinations. The primers are at 10 $\mu$ M each and are mixed in equal parts to create the 10 $\mu$ M combinations.

We use the bacteria/archaeal primers 515F/806R with an inhouse barcode system designed by Aaron Darling.

More about the primers can be found here:

Caporaso JG, Lauber CL, Walters WA, Berg-Lyons D, Huntley J, Fierer N, Owens SM, Betley J, Fraser L, Bauer M, Gormley N, Gilbert JA, Smith G, Knight R. 2012. Ultra-high-throughput microbial community analysis on the Illumina HiSeq and MiSeq platforms. ISME J.

---

### 1X Recipe:

22.5  $\mu$ L of Invitrogen Platinum SuperMix

1.25  $\mu$ L of Primer Mix (10 $\mu$ M)

0.62  $\mu$ L of PNA chloroplast blocker ( $\sim$  50 $\mu$ M)

0.62  $\mu$ L of PNA mitochondria blocker ( $\sim$  50 $\mu$ M)

1.0-5.0  $\mu$ L of template DNA

**Total:** 26-30 $\mu$ L

---

### Protocol:

#### *PCR*

1. Amplify samples in triplicate, meaning each sample will be amplified in 3 replicate 25  $\mu$ L PCR reactions.
2. Combine the triplicate PCR reactions for each sample into a single volume. Combination will result in a total of 75  $\mu$ L of amplicon for each sample. Do NOT combine amplicons from different samples at this point.
3. Run 5  $\mu$ L of amplicons for each sample on an agarose gel to verify amplification. Expected band size for roughly 300-400 bp

#### *Clean-up*

4. Clean PCR reactions with magnetic beads. See "Magnetic-bead PCR Cleaning Protocol" for an example protocol.
5. Quantify amplicons with Qubit.

#### *Pooling*

6. Combine an equal amount of amplicon from each sample into a single, sterile tube. Generally 240 ng of DNA per sample are pooled. However, higher amounts can be used if the final pool will be gel-isolated or when working with low biomass samples. *Note:* When working with multiple plates of samples, it is typical to produce a single tube of amplicons for each plate of samples.
7. If spurious bands were present on gel (in step 3), half of the final pool can be run on a gel and then gel-extracted to select only the target bands.\*

\*We routinely see high molecular-weight bands or smears in the gel. We use the Pippin prep gel from Sage Science (<http://www.sagescience.com/products/pippin-prep/>) to select the correct band size.

8. Measure concentration and 260/280 of final pool that has been cleaned. For best results the 260/280 should be between 1.8-2.0.
9. Send an aliquot for sequencing along with sequencing primers.

---

**Thermocycler Conditions:**

*For 96 well thermocyclers:*

1. 94°C 3 minutes
2. 94°C 45 seconds
3. 50°C 60 seconds
4. 72°C 90 seconds
5. Repeat steps 2-4 35 times
6. 72°C 10 minutes
7. 4°C HOLD

*For 384 well thermocyclers:*

1. 94°C 3 minutes
2. 94°C 60 seconds
3. 50°C 60 seconds
4. 72°C 105 seconds
5. Repeat steps 2-4 35 times
6. 72°C 10 minutes
7. 4°C HOLD